


Nebulizer Piston Type


Model : NE100



- MMAD <4 μm ; Fine Particle Dose (FPD): 75-80%
- Consistent fine particle size for efficient respiratory treatment
- Powerful piston compressor
- Mouthpiece and masks for adult and child included
- Economic design



Consistent fine particle size for efficient respiratory treatment



Adult and child masks




Nebulizer kit (N2)



Aerosol Performance Tested with 0.9% saline			
Cascade Impactor		Malvern Spraytec	
MMAD	FPD	Dv(50)	FPD
3.01 μm	75.02%	4.71 μm	53.46%

FPD: Fine particle dose, the percentage of particle size less than 5.0 μm



Model	Qty per carton
NE100	-----

● Storage and Transportation Condition
 -20°C to 60°C (-4°F to 140°F)
 10% - 90% RH, 700 - 1060hPa